

ZOLL® AED Plus®

Automated External Defibrillator

Real Help. Real Fast.

The ZOLL AED Plus® is designed to help rescue more victims of sudden cardiac arrest. Easy to use, highly durable, and packed with life-saving technology.



Clear Voice Guidance
Step-by-step voice prompts lead you through the rescue.



Real CPR Help®
Provides real-time feedback on CPR depth and rate.



Long Battery Life
Up to 5 years of standby life.



Rugged & Reliable
Built to perform when and where it matters most.



Easy to use. Hard to ignore.
Ready when every second counts.



PRODUCT VIEWS



Front View



Side View



Back View



Bottom View

KEY FEATURES



One-Piece Electrode Pad
Electrode pad with Real CPR Help® for faster, easier placement.



Simple, Clear Prompts
Calm voice instructions guide you every step of the way.



Automatic Self-Test
Performs daily self-checks to ensure it's always ready.



Child & Adult Use
Pediatric capability with the push of a button.



Durable Design
Built to withstand harsh environments (IP55 rating).

SPECIFICATIONS

Dimensions	13.3 cm (H) x 24.1 cm (W) x 29.2 cm (D)
Weight	3.1 kg (with battery and pads)
Battery Type	10.8 V lithium manganese dioxide (lithium primary)
Battery Life	Up to 5 years (standby) or 225 shocks / 13 hours of operation
Electrode Pad Shelf Life	Up to 5 years
Operating Temperature	0°C to 50°C
Ingress Protection	IP55
Warranty	7 years

USAGE STEPS



Turn On
Press the ON/OFF button to start.



Apply Pads
Place the electrode pad on the patient's bare chest.



Follow Prompts
Listen carefully to the clear voice instructions.



Shock If Advised
If a shock is needed, press the flashing shock button.



Continue CPR
Continue CPR as instructed until help arrives.



Save More Lives
Advanced technology that helps rescuers perform with confidence.



Real CPR Help®
Provides real-time feedback on CPR quality to help improve outcomes.



Always Ready
Daily self-tests ensure your AED is ready when you need it.

ZOLL®

Advancing Resuscitation.
Today. Tomorrow. Together.

ZOLL AED PLUS

Automated External Defibrillator

The ZOLL AED Plus is a trusted, easy-to-use AED designed for rapid response in sudden cardiac arrest (SCA) emergencies. It combines real-time CPR feedback, intuitive voice prompts, and advanced technology to help rescuers deliver high-quality care with confidence.



Includes

- ✓ ZOLL AED Plus Unit
- ✓ CPR-D-padz® Adult/Child Electrode Pads
- ✓ Duracell Lithium Battery (10-year shelf life)
- ✓ Soft Carry Case
- ✓ User Guide
- ✓ Ready-to-use Packaging



Real CPR Help®

Provides real-time feedback on chest compression rate and depth.



Clear Voice Prompts

Step-by-step audio instructions guide users through the rescue process.



Easy to Use

One-button operation and intuitive design make it ideal for any responder.



Durable & Reliable

Built to withstand harsh environments (IP55 rating for dust and water resistance).



Proven Technology

Trusted by professionals and lay rescuers worldwide.



Long Battery Life

5-year battery life for readiness when it matters most.



For Professionals
and the Public



Trusted by
Millions



Designed for
Real Emergencies



Helping Save
More Lives