

MEDICAL EXAMINATION BED

Professional Clinical Examination Couch

MEDICAL EXAMINATION BED



6 cm Thick Sponge with PVC Cover



Adjustable Backrest Up to 60°



Sturdy Metal Frame with Anti-Rust Coating



Reinforced Structure for Extra Stability



Non-Slip Rubber Foot Caps



COMFORTABLE MEDICAL EXAMINATION BED

This bed is specially designed to provide maximum comfort and stability during the medical examination.



ADJUSTABILITY

The backrest tilt can be adjusted up to 60 degrees; To make it easier for the doctor to conduct the examination easily.



DURABLE METAL FRAME

The metal frame coated with an anti-rust material ensures durability and endurance.



COMFORTABLE & CLEAN SURFACE

The 6cm thick sponge surface is wrapped with anti-bacterial PVC, ensuring comfort and cleanliness.



IDEAL DIMENSIONS

With dimensions of 185*62*72cm, the table can fit into a variety of spaces.



WEIGHT-BEARING CAPACITY

This table can bear a weight of up to 160 kilograms, making it able to accommodate a variety of individuals.

Product Overview

This medical examination bed is specially designed to provide maximum comfort, stability, and durability for hospitals, clinics, and healthcare facilities.

Key Features

- Comfortable medical examination bed
- Adjustable backrest up to 60°
- Durable anti-rust coated metal frame
- 6 cm thick sponge surface with anti-bacterial PVC
- Ideal dimensions: 185 × 62 × 72 cm
- Weight-bearing capacity up to 160 kg

Technical Specifications

Specification	Details
Overall Size	185 × 62 × 72 cm

Backrest Adjustment	0° – 60°
Frame	Anti-Rust Coated Metal
Surface	Anti-Bacterial PVC
Foam Thickness	6 cm
Weight Capacity	Up to 160 kg

ENTER Medical Furniture Solutions

