


YA-ET15F – Examination Couch with Back Support

Elegant • Heavy-Duty • Clinical Grade



Product Overview

The YA-ET15F Examination Couch is a durable and elegant medical examination table designed for hospitals, clinics, diagnostic centers, and healthcare facilities. Built with a strong powder-coated steel structure and premium upholstered cushioning, it provides excellent patient comfort, stability, and long-lasting performance.

Key Features

- Heavy-duty square tube steel construction
- Epoxy powder-coated metal frame
- Adjustable backrest with ratchet mechanism
- High-density foam upholstered top
- Easy-to-clean synthetic leather upholstery
- Rubber foot caps for stability
- Knock-down structure for transport and storage
- Optional paper roll holder available

Technical Specifications

Specification	Details
Model	YA-ET15F
Product Type	Examination Couch
Overall Size	180 × 60 × 70 cm
Backrest Adjustment	0° – 65°
Frame Material	Mild Steel Square Tube
Finish	Epoxy Powder Coated
Upholstery	High Quality Synthetic Leather
Weight Capacity	Up to 135 kg

Applications

Suitable for hospitals, medical clinics, examination rooms, physiotherapy centers, healthcare facilities, and diagnostic centers.

Professional Medical Furniture Solution